ThermaCare HeatWraps, Multi-Purpose

Muscle Pain Therapy

This is the most current labeling information and may differ from labels on product packaging. If there are any differences between this website labeling and product packaging labeling, this website labeling should be regarded as the most current.

Uses
 Provides heat therapy for temporary relief of minor muscular and joint aches and
pains associated with overexertion, strains, sprains, and arthritis.
 Warnings THIS PRODUCT CAN CAUSE BURNS 55 OR OLDER: YOUR RISK OF BURNING INCREASES AS YOU AGE. WHEN APPLIED OVER THE LOWER BACK/HIP, WEAR THERMACARE OVER A LAYER OF CLOTHING, NOT DIRECTLY AGAINST YOUR SKIN. WHEN APPLIED OVER THE KNEE OR ELBOW, WEAR THERMACARE OVER A TOWEL OR CLOTH SUCH AS A WASHCLOTH, NOT DIRECTLY AGAINST YOUR SKIN. DO NOT WEAR WHILE SLEEPING. CHECK SKIN FREQUENTLY DURING USE. IF YOU FIND IRRITATION OR A BURN, REMOVE PRODUCT IMMEDIATELY.
 Do not use on the face/head, breasts, lower abdomen, fingers, groin, feet or toes on the back of the knee or inside the bend of the arm or underarm if heat cell contents leak and/or wrap is damaged or torn for more than 8 hours in a 24-hour period with pain rubs, medicated lotions, creams or ointments on unhealthy, damaged or broken skin on areas of bruising or swelling that have occurred within 48 hours on areas of the body where you can't feel heat with other forms of heat on people unable to remove the product, including children, infants, and some elderly on people unable to follow all use instructions
Ask a doctor before use if you have DIABETES poor circulation or heart disease rheumatoid arthritis or are pregnant
 When using this product check skin frequently for signs of burns or blisters — if found, stop use if product feels too hot — stop use or wear over clothing do not place extra pressure over the product such as leaning against hard surfaces, or wearing under tight clothing, a tight waistband or belt to avoid increased risk of burns, remove product before undergoing any MRI (magnetic resonance imaging) procedures

ThermaCare HeatWraps, Multi-Purpose Muscle Pain Therapy

This is the most current labeling information and may differ from labels on product packaging. If there are any differences between this website labeling and product packaging labeling, this website labeling should be regarded as the most current.

Stop use and ask a doctor

- if pain worsens or persists for more than 7 days
- if you experience any discomfort, burning, swelling, rash or other changes in your skin that persist where the wrap is worn

Keep out of reach of children and pets. Each heat cell contains iron (~2 grams):

- Can be harmful if ingested. If ingested, rinse mouth with water and call a Poison Control Center immediately.
- If heat cell contents come in contact with your skin or eyes, remove wrap, rinse affected area with water, and seek professional assistance immediately.

DO NOT MICROWAVE or attempt to reheat as wrap could catch fire

Directions

- Tear open pouch when ready to use (it may take up to 30 minutes to reach therapeutic temperature).
- Peel away paper to reveal adhesive side.
- Place on pain area with adhesive against the skin. Refer to DO NOT USE section for areas where not to apply.
- Attach firmly.
- Be careful when applying to the wrist <u>do not</u> overlap the heat cells.
- Place wrap flat. Do not fold, overlap or bunch the heat cells as your risk of burn is increased.
- If 55 and older:
 - When applied over the lower back/hip, wear ThermaCare over a layer of clothing, not directly against your skin.
 - When applied over the knee or elbow, wear ThermaCare over a towel or cloth such as a washcloth, not directly against your skin.
- Wear up to 8 hours. <u>Do not</u> wear for more than 8 hours in any 24-hour period.
- Discard after use.

Questions?

1-800-323-3383